



Bottle, sampling, presterilized, round, 50/pk

Продукт №: 2495050
BGN Цена: Свържете се с нас
Не е показана дата на доставка

Polystyrene, disposable, transparent, non-fluorescent, break-resistant, presterilized bottles with caps. Bottles have 100 mL fill-to line; capacity of 175 mL. Bottles used with medium containing MUG or other fluorescing reagent. Work well with P-A Broth Ampules and PathoScreen Medium.

Спецификации

Dechlorinating Agent?:	No
Обем / опаковка, размер:	50 Stück