



Bottle Carrier, Linear

Продукт №: 2094300
BGN Цена: Свържете се с нас
Не е показана дата на доставка

Linear, low-density polyethylene. Securely holds up to six standard 0.5 L glass bottles, pint-size reagent bottles, and most other 500 mL glass or plastic bottles—making breakage almost impossible. Lightweight, chemical-resistant, and easy to clean. Dimensions: 34.9 x 24.8 x 27.6 cm (14 x 10 x 11"). Maximum bottle size it will hold: 73 mm in diameter, 178 mm in height.

Спецификации

Дълбочина:	27.6 cm
Дължина:	349 mm
Обем / опаковка, размер:	1 each
Ширина:	10 Inches