




Filter, acrodisc, pore size 0.2 µm

Продукт №:	2415600
BGN Цена:	Свържете се с нас
Не е показана дата на доставка	

For filtration of dilute biological fluids. Sterilized by gamma irradiation to eliminate potential contamination by EtO residuals. Inlet/Outlet connections: 13/25 mm: Female Luer-Lok inlet, standard male slip luer outlet. Pore size: 0.2 µm. Diameter: 25mm.

Спецификации

Filter Type:	Filter, acrodisc
Packing:	NA
Диаметър:	NA
Материал:	Polysulfone
Пори, размер:	0.2 µm