



Filter, membrane, black, pore size 0.45 µm, diameter 47 mm, 200/pk, w/out pads

Продукт №:

2490300

BGN Цена:

Свържете се с нас

Не е показана дата на доставка

Provides contrast for microbiological and particulate analysis. Modified polyethersulfone creates dark background for excellent contrast useful in enumerating opaque colonies and light-covered particulate. White grid line aids enumeration. Certified for membrane filtration analysis. Sterile and individually packaged. Pore size 0.45 µm. Diameter: 47 mm.

Спецификации

Filter Type:	Filter
Диаметър:	47 mm
Материал:	Modified polyethersulfone
Обем / опаковка, размер:	200 /pk
Пори, размер:	0.45 µm