



## Pitcher, graduated, 1000 cc

Продукт №: 2612853  
BGN Цена: Свържете се с нас  
Не е показана дата на доставка

Chemical-resistant, shatter-proof, tall-form pitchers with well-designed spouts, and large handles. Metric graduations are molded onto outside wall of each pitcher. Polypropylene, steam autoclavable at 250°F (121°C).

---

### Спецификации

Graduations:	50 mL
Дълбочина:	167 mm
Капацитет:	1000 mL
Обем / опаковка, размер:	1 each