



Pitcher, graduated, 2000 cc

Продукт №: 2612854
BGN Цена: Свържете се с нас
Не е показана дата на доставка

Chemical-resistant, shatter-proof, tall-form pitchers with well-designed spouts, and large handles. Metric graduations are molded onto outside wall of each pitcher. Polypropylene, steam autoclavable at 250°F (121°C).

Спецификации

Graduations:	100 mL
Дълбочина:	205 mm
Капацитет:	2000 mL
Обем / опаковка, размер:	1 each