



Filter, membrane, pore size 0.2 μm , diameter 47 mm, 200/pk (Gelman)

Продукт №:

2618800

BGN Цена:

Свържете се с нас

Не е показана дата на доставка

For isolation and enumeration of stressed organisms (*Pseudomonas* species). High retention and flow rate. Certified for use in membrane filtration analysis. White polyethersulfone media. STerile and individually packaged. Pore size: 0.2 μm . Diameter: 47mm.

Спецификации

Filter Type:	Filter
Диаметър:	47 mm
Материал:	Polyethersulfone
Обем / опаковка, размер:	200 /pk
Пори, размер:	0.2 μm